

# GCare *Lipid*

Point of Care Lipid Profile



Customer Care




**Total Kolesterol**  
**Trigliserid**  
**HDL-C**  
**LDL-C (Kalkulasi)**  
**Glukosa**



 **GCMS**

## GCare *Glucose*

**GCare Glucose Test Strip**



**STRIP BESAR,  
MUDAH DIPEGANG,  
NYAMAN, HIGIENIS**

Pengujian harus selesai dalam waktu 1 menit setelah memasukkan tes strip

GCM-110C

Mudah  
Digunakan




Printer Eksternal



Self Optical Check

# GCare *Lipid*



## Blood glucose test strip

|                          |                                                          |
|--------------------------|----------------------------------------------------------|
| <b>Nama Produk</b>       | GCare Glucose Test Strip                                 |
| <b>Packing unit</b>      | 50 tes / BOX                                             |
| <b>Model</b>             | GBP-110                                                  |
| <b>Suhu Penyimpanan</b>  | 1~32°C (34°F~90°F)                                       |
| <b>Material Provided</b> | GCare Glucose Test Strip x<br>50 Petunjuk penggunaan x 1 |

## GCare Lipid analyzer

|                           |                                                                                         |
|---------------------------|-----------------------------------------------------------------------------------------|
| <b>Nama Produk</b>        | GCare Lipid                                                                             |
| <b>Memori</b>             | Menyimpan 1000 data Kolesterol<br>Menyimpan 500 data Glukosa. GCM-                      |
| <b>Model</b>              | 110, GCM-110W, GCM-110B                                                                 |
| <b>Berat</b>              | 120g                                                                                    |
| <b>Ukuran Sampel</b>      | Kolesterol 40µl / Glukosa 0.5                                                           |
| <b>Dimensi</b>            | 130 x 65 x 30 mm <sup>3</sup>                                                           |
| <b>Hematocrit</b>         | Kolesterol 30~50% / Glukosa 20~60%                                                      |
| <b>Waktu</b>              | Kolesterol 150 detik / Glukosa 5 detik                                                  |
| <b>Daya</b>               | (AAA) 1.5V X 4                                                                          |
| <b>Kisaran Pengukuran</b> | TC 100 ~ 400 mg/dL<br>TG 50 ~ 600 mg/dL<br>HDL 25 ~ 100 mg/dL<br>Glukosa 10 ~ 900 mg/dL |
| <b>Suhu</b>               | Kolesterol 18~30° Glukosa 4~40°                                                         |
| <b>Unit</b>               | mg/dL or mmol/L                                                                         |
| <b>Penyimpanan</b>        | -10~60° (14°~140°)                                                                      |
| <b>Kelembapan</b>         | 10~75% R.H.                                                                             |

## Lipid profile test strip

|                          |                                                                                                                 |
|--------------------------|-----------------------------------------------------------------------------------------------------------------|
| <b>Nama Produk</b>       | GCare Lipid Profile Test Strip                                                                                  |
| <b>Packing unit</b>      | 10 tes / BOX                                                                                                    |
| <b>Model</b>             | GLP-110                                                                                                         |
| <b>Suhu Penyimpanan</b>  | 2~30°C (36°F~86°F)                                                                                              |
| <b>Material Provided</b> | GCare Lipid Profile Test Strip x<br>10 Alat tetes volume tetap x<br>10 Petunjuk penggunaan x 1<br>Kode chip x 1 |



Diimpor Oleh :

Diproduksi Oleh :



**PT. GANDASARI EKASATYA**  
Medical, Laboratory, Scientific Supply & General Trading



**GC**  
Medical Science

Jl. Kemanggisan Ilir Raya Blok C NO. 15B  
Anggrek Neli Murni Jakarta 11480  
Telp : (+62-21) 8400484, Hotline : +62 811 9940 129  
E-mail : support@gandasari-ekasatya.co.id

(16924) 107, Ihyeon-ro 30beon-gil,  
Giheung-gu, Yongin-si, Gyeonggi-do, Korea  
TEL +82-31-260-0965  
www.greencrossms.com